

store.independenceaustralia.com
customerservice@iagroup.org.au
T 1300 788 855 F 1300 788 811

independence
AUSTRALIA
Your way.

Compression Products Stockings & Accessories

2021/22

Shop for a
good cause 





Venous leg ulcer therapy

Advanced woundcare and compression therapy

With products - **Cutimed® Sorbion®** and **JOBST® FarrowWrap®**, JOBST® enables effective management of venous leg ulcers whilst offering ease of use for clinicians and patients.

The perfect combination for optimal wound healing and patient empowerment.

For more information please contact Independence Australia on 1300 788 855 or visit store.independenceaustralia.com



Contents

Our Services	4
How To Place An Online Order	5
Understanding Compression Stockings	6
Sizing Charts	7
How to Measure for Compression Stockings	10
How to Apply & Remove Compression Stockings	12
Anti Bacterial Stockings	16
Anti Embolism Stockings	18
Body Compression Garments	22
Compression Socks	24
Compression Stockings	28
Lymphoedema Treatments	44
Stocking Applicators & Aids	46
Ulcer Treatments	48

our Services

Independence Australia specialises in delivering quality health and wellness products to people in their homes to help them enjoy a better quality of life.

We stock over 12,000 products across all the major brands and offer expert advice to ensure you have the best products to meet your needs.

Why choose us?

Competitive prices

- Stretch your dollar further
- Exclusive product promotions to guarantee value for money

Prompt and professional delivery

- Warehouses in each state for quick delivery nationwide
- Free delivery for all online orders over \$50

Easy payment options

- A 14-day payment plan
- A variety of payment methods are available

Professional service for your convenience

- Customer Service is available Monday to Friday 8am to 7pm with extended hours during daylight savings
- A team of professional staff to provide expert advice and support



Quick and reliable delivery



Discreet packaging available



Free delivery for online orders over \$50

Conditions apply.
Visit store.independenceaustralia.com

Shop for a good cause

As a social enterprise, the income we make is invested back into the vital services we provide to help people with disabilities. *Your purchases help us to enhance people's lives.*

Order today

 store.independenceaustralia.com

 1300 788 855

 customerservice@iagroup.org.au

 1300 788 811

Payment methods

Credit Card



BPay



Post billpay



PayPal
(Online orders)



Independence
AUSTRALIA
Your way.

How to place an online order with Independence Australia


If you're new to online shopping, here's how to shop with Independence Australia in a few simple steps.

1



Create Shopping Basket


To get started, type store.independenceaustralia.com into your web browser.

1. Use **search bar** at the top of the screen 

2. **Select item** with [Add to Basket](#)

3. **Enter quantity** PKT Qty

4. **Click** [Add to Basket](#) to move item to shopping basket. Repeat steps until shopping is complete.

5. **Click**  to review your shopping list.

6. **Click** [Checkout Now](#) when product selection is complete.

2



Confirm Shopping

7. **Review** order summary.

8. **Add your personal details**, choose a delivery option and confirm any delivery instructions.

9. **Click** [Next](#)

3



Checkout

10. **Choose a payment method** and enter your credit card or PayPal details.

11. **Tick** ☐ I accept the Terms and Conditions*

12. **Click** [Place Order](#) to complete your transaction.

13. **Look out for an email** confirming your order details.


Tips for Ordering Online

- You'll need an **email address** to place an order
- Choose **"Discreet Packaging"** to keep your delivery box details private
- Hit **"Back"** button at any time to go back a step
- **Set up an account** so you can reorder products easily and track order status

Fact Sheet For Understanding Compression Stockings

What are compression stockings?

Snug-fitting, stretchy socks that gently squeeze your leg with graduated compression or pressure. They will be tighter around your ankle and get looser as they sit up higher on your leg.



Who should wear compression stockings?

Recommended for anyone who has venous disease (abnormal circulatory condition), Lymphoedema, spends time sitting, standing or in a sedentary position including people who have difficulty moving their legs.

Compression stockings should be worn with caution for people with:

- skin infections
- weeping wounds
- allergies to garment fabric
- impaired sensitivity of the limb
- immobility (confined to bed)

Why wear compression stockings?

Compression stockings improve blood flow by helping to prevent pain, swelling or build-up of fluid in your legs. They can also lower the chances of deep vein thrombosis (DVT) and other circulation problems.

What is mmHG?

mmHG stands for millimetres of mercury and is the level of pressure or compression of the stocking, usually expressed as a number range.

There are 4 main compression levels:

- 8-15mmHG [Mild Compression]
- 15-20mmHG [Moderate Compression]
- 20-30mmHG [Firm Compression]
- 30-40mmHG [Extra Firm Compression]

The higher the numbers, the higher the compression. The amount of pressure at the top of the stocking is the lower number in the mmHG range, while the higher number is the amount of pressure at the foot of the stocking. Stockings with higher levels of compression (20-30mmHG to 30-40mmHG) will need a prescription from a doctor.

Where can you buy stockings from?

Independence Australia stocks hundreds of options for compression stockings online at store.independenceaustralia.com.

Tips & Advice

- Stockings should feel snug, but not painfully tight. They need to be firm fitting to do their job and so they may feel tight at first.
- It is normal for stockings to be pulled up or repositioned a few times a day.
- Open-toe stockings are often easier to get on and off but are otherwise the same.
- Stockings need to be replaced every 3 to 6 months, or when they start to sag or become loose.
- Compression stockings should be worn even during hot summer months, when legs need the support more than ever.
- Handwash gently in warm water with a mild soap or detergent, or in a clothes bag on a gentle machine cycle. Do not use bleach, fabric softeners or hot water.
- Before reordering replacements, remeasure legs and if measurements are greater or less than 3cm from original measurement, you may need a new size.

To access a wide range of compression stockings, you can visit store.independenceaustralia.com for quick and reliable home delivery

© Independence Australia Group 2021. No part of this information may be reproduced in whole or in part in any manner without permission.

Compression Stocking Sizing Charts

For your easy reference, the following are the sizing charts for our most popular compression stocking brands. Please refer to a manufacturer's website or packaging for the most up to date sizing information.

JOBST				
SIZE	ANKLE	CALF	THIGH	HIP
JOBST ULCERCARE				
SMALL	18cm - 21cm	29cm - 35cm		
MEDIUM	21cm - 24cm	34cm - 41cm		
LARGE	24cm - 27cm	40cm - 47cm		
EXTRA LARGE	27cm - 30cm	46cm - 54cm		
XX LARGE	30cm - 33cm	48cm - 56cm		
XXX LARGE	33cm - 36cm	50cm - 58cm		
XXXX LARGE	36cm - 39cm	52cm - 60cm		
JOBST ULTRASHEER, FOR MEN, FOR MEN CASUAL, ACTIVEWEAR, SPORT AND RELIEF (MATERNITY)				
SMALL	18 - 21cm		40cm - 62cm	81cm - 132cm
MEDIUM	21 - 25cm		46cm - 70cm	86cm - 142cm
LARGE	25 - 29cm		54cm - 78cm	91cm - 152cm
EXTRA LARGE	29 - 33cm		60cm - 81cm	112cm - 180cm
JOBST ULTRASHEER, FOR MEN, FOR MEN CASUAL, ACTIVEWEAR, SPORT AND RELIEF (UNISEX)				
SMALL	18cm - 21cm	28cm - 38cm	40cm - 62cm	71cm - 117cm
MEDIUM	21cm - 25cm	30cm - 42cm	46cm - 70cm	76cm - 127cm
LARGE	25cm - 29cm	32cm - 46cm	54cm - 78cm	81cm - 137cm
LARGE FC	25cm - 30cm	46cm - 61cm		
EXTRA LARGE	29cm - 33cm	34cm - 50cm	60cm - 81cm	102cm - 166cm
EXTRA LARGE FC	30cm - 36cm	46cm - 61cm		

*FC = Full Calf

Compression Stocking Sizing Charts

VENOSAN					
SIZE		ANKLE	CALF	SHOE SIZE	
KNEE HIGH SUPPORT SOCKS (SILVERLINE & MICROFIBRE) CLASS I 15-20 mmHg CLASS II 20-30mmHg					
MALE	SMALL	19cm -21cm	28cm - 38cm	7 & Under	
	MEDIUM	21.5cm -24cm	33cm - 43.5cm	7 1/2 - 10	
	LARGE	24cm - 26.5cm	35.5cm - 46cm	10 1/2 - 12	
	EXTRA LARGE	26.5cm - 32cm	38.5cm - 51.5cm	12 1/2 +	
	XXL - SILVERLINE ONLY	32cm - 36cm	41cm - 53cm	13-15	
FEMALE	SMALL		26.5cm - 37cm	4-5	
	MEDIUM		29cm - 39.5cm	5 1/2 - 7 1/2	
	LARGE		32cm - 43cm	8-10 1/2	
	EXTRA LARGE		34.5cm - 47cm	11-13	
	XXL - SILVERLINE ONLY		37cm - 49.5cm	14-15	
SIZE		ANKLE	BELOW KNEE	MID CALF	THIGH HIGH
VENOSAN STANDARD SIZING CHART FOR 4000/5000/7000 CLASS 1 18-21MMHG CLASS II 23-32MMHG CLASS III 34-46MMHG					
SMALL		20cm - 23cm	28cm - 38cm	29cm - 39cm	48cm - 63cm
MEDIUM		23cm - 26cm	32cm - 41cm	33cm - 42cm	54cm - 67cm
LARGE		26cm - 29cm	35cm - 44cm	35cm - 45cm	58cm - 72cm
EXTRA LARGE		29cm - 31cm	38cm - 46cm	39cm - 47cm	63cm - 74cm
LENGTH		ID = Below Knee	IK = Pantyhose	IG = Thigh High	
SHORT		35cm - 40cm	67cm - 74cm	63cm - 72cm	
REGULAR/LONG ALL SERIES		40cm - 44cm	74cm - 85cm	68cm - 78cm	
*Available in Custom Made					
VENOSAN STANDARD SIZING CHART FOR 6000 CLASS 1 18-21MMHG CLASS II 23-32MMHG CLASS III 34-46MMHG					
SMALL		20cm - 22cm	29cm - 38cm	30cm - 39cm	49cm - 63cm
MEDIUM		23cm - 25cm	32cm - 40cm	33cm - 42cm	54cm - 66cm
LARGE		26cm - 28cm	35cm - 43cm	36cm - 44cm	58cm - 70cm
EXTRA LARGE		29cm - 31cm	38cm - 46cm	39cm - 47cm	63cm - 74cm
LENGTH		ID = Below Knee	IK = Pantyhose	IG = Thigh High	
SHORT		34cm - 39cm	67cm - 74cm	59cm - 67cm	
REGULAR/LONG ALL SERIES		39cm - 44cm	74cm - 85cm	68cm - 78cm	

HEALTH SUPPORT SYSTEM

SIZE	ANKLE	CALF	CALF HEIGHT (KNEE TO HEEL)		LENGTH
			BEIGE HOISERY	BLACK SOCK	
A	18cm -20cm	25.5cm - 33cm	To 38cm Over 38cm	To 40.5cm Over 40.5cm	S R
B	20cm -22.5cm	30.5cm - 38cm	To 40.5cm Over 40.5cm	To 43cm Over 43cm	S R
C	22.5cm -25.5cm	35.5cm - 43cm	To 43cm Over 43cm	To 45.5cm Over 45.5cm	S R
D	25.5cm -28cm	40.5cm - 48.5cm	To 45.5cm Over 45.5cm		S R
E	28cm -30.5cm	45.5cm - 53cm	To 45.5cm Over 45.5cm		S R
F	30.5cm - 33cm	50cm - 58.5cm	To 45.5cm Over 45.5cm		S R
G	33cm - 35.5cm	55.5cm - 66cm	To 45.5cm Over 45.5cm		S R

COMPRINET

SIZE	LENGTH	WIDTH
PRO THIGH HIGH		
SMALL	Regular 64cm - 80cm	23cm - 30cm
SMALL	Long 80cm - 97cm	23cm - 30cm
MEDIUM	Regular 64cm - 80cm	30cm - 38cm
MEDIUM	Long 80cm - 97cm	30cm - 38cm
LARGE	Regular 64cm - 80cm	38cm - 46cm
LARGE	Long 80cm - 97cm	38cm - 46cm
EXTRA LARGE	Regular 64cm - 80cm	46cm - 53cm
EXTRA LARGE	Long 80cm - 97cm	46cm - 53cm
PRO KNEE HIGH		
SMALL		23cm - 30cm
MEDIUM		30cm - 38cm
LARGE		38cm - 46cm
EXTRA LARGE		46cm - 53cm

THROMBEXIN

SIZE	ANKLE	CALF	COLOUR CODE
KNEE LENGTH STOCKINGS			
EXTRA SMALL	17cm - 19cm	26cm - 32cm	
SMALL	20cm - 22cm	30cm - 38cm	
MEDIUM	23cm - 25cm	34cm - 44cm	
LARGE	26cm - 28cm	38cm - 48cm	
LARGE / EXTRA WIDE		44cm - 55cm	
EXTRA LARGE	29cm - 31cm	44cm - 55cm	
THIGH LENGTH STOCKINGS			
EXTRA SMALL	17cm - 19cm	40cm - 50cm	
SMALL	20cm - 22cm	48cm - 60cm	
SMALL / EXTRA WIDE		56cm - 70cm	
MEDIUM	23cm - 25cm	56cm - 70cm	
MEDIUM / EXTRA WIDE		64cm - 80cm	
LARGE	26cm - 28cm	64cm - 80cm	
LARGE / EXTRA WIDE		72cm - 90cm	
EXTRA LARGE		72cm - 90cm	
EXTRA LARGE / EXTRA WIDE	29cm - 31cm	80cm - 100cm	
THIGH LENGTH STOCKINGS WITH WAIST ATTACHMENT			
SMALL	20cm - 22cm	48cm - 60cm	
MEDIUM	23cm - 25cm	56cm - 70cm	
LARGE	26cm - 28cm	64cm - 80cm	
LARGE / EXTRA WIDE		72cm - 90cm	



Fact Sheet For How to Measure for Compression Stockings

Identifying the correct size of compression stockings is essential to ensure you have a good fit, feel comfortable and receive the optimal benefits. If you have never worn compression stockings before, or are unsure you have the correct size, this guide will help you find the right fit.

Getting started

You will need:

- Tape measure (preferably soft tape or cloth)
- Sizing chart for your chosen stocking brand

Remember

- There is no universal standard for stocking sizes. Different brands use different sizing charts. For example, a medium size in one brand may not be a medium in another brand.
- This is not a guide for travel or flight socks which are class 1 stockings (lowest in compression to prevent DVT) and require only basic information such as height, weight or foot size for measurement.






Tips


- Measure earlier in the day before swelling builds in the legs.
- Measure next to bare skin and without shoes on.
- Always place foot flat against the floor and the leg at a right angle to the foot.
- Take length measurements from the ground up.
- Always pull the tape snugly and without constriction around a leg.
- If a tape measure is unavailable, use a piece of string and measure individual lengths with a ruler.
- If you are in between two sizes, choose the size that suits the ankle circumference best.
- If the foot size of a stocking is either too large or too small, try a stocking with an open toe.

Determine measurement points

To start, take note of which area of the leg needs to be measured for your preferred stocking style. Compression stockings require you to measure the circumference at specific points on the leg and the length of the leg.

Stocking style	Area of Leg				
	Ankle circumference	Calf circumference	Calf length	Thigh circumference	Leg length
 Knee high	✓	✓	✓		
 Thigh high	✓	✓		✓	✓
 Waist high	✓	✓		✓	✓

Start Measuring




Ankle measurement is _____

Measure your ankle circumference

- Place a measuring tape at the narrowest part of the ankle, just above the ankle bone, and measure the circumference.
- The ankle must always be your first measurement and is the most important, as the graduated compression begins at the ankle.



Ankle measurements may differ from leg to leg. In this instance note the larger size.




Calf circumference measurement is _____

Measure your calf circumference

- Place a measuring tape at the largest part of the calf and measure the circumference.



Choose a size that puts your calf measurement in the middle of the sizing range to avoid stockings that are too tight and make the top band bind or roll, or too large making the stocking slip or fall.




Calf length measurement is _____

Measure your calf length

- Sit on a chair and measure the distance from the floor to just below the bend in your knee.




Do not wear shoes.



Thigh circumference measurement is _____

Measure your thigh circumference

- Find the widest part of your thigh, directly under your buttocks and measure the circumference.



Leg length measurement is _____

Measure your leg length

- Measure the distance from the floor behind the heel to the bottom of your buttocks (top of thigh).



Do not wear shoes. Keep leg straight.

Final Step

You can now compare your measurements with the sizing chart of your preferred stocking brand to determine the best size and length for your compression stockings.

How to Use a Sizing Chart

As an example, to measure for a knee-high stocking, the measurement points are ankle circumference, calf circumference and calf length. If measuring results in **20cm** for ankle circumference, **37cm** for calf circumference and **39cm** for calf length, your size is determined as **A**.

Sample Sizing Chart

Ankle circumference	Small 18 - 21cm		Medium 21 - 25cm		Large 25 - 30cm	
Calf circumference	Small 28 - 39cm		Medium 31 - 45cm		Large 35 - 51cm	
Calf length	Short <40cm	Long >40cm	Short <40cm	Long >40cm	Short <40cm	Long >40cm
Thigh circumference	Small 45 - 65cm		Medium 48 - 75cm		Large 52 - 85cm	
Leg length	Short <70cm	Long >70cm	Short <70cm	Long >70cm	Short <70cm	Long >70cm
Brand size	A	B	C	D	E	F

Note: The sizing chart of your body

Note: The sizing chart of your stocking brand will look different to this table, but they all follow the same principles. Refer to a brand's website or packaging to find a sizing chart.

To access a wide range of compression stockings visit

store.independenceaustralia.com for quick and reliable home delivery

© Independence Australia Group 2021. No part of this information may be reproduced in whole or in part in any manner without permission.

Fact Sheet For How to Apply & Remove Compression Stockings

Compression stockings are often difficult to apply and remove as the compression is greatest at the ankle, requiring the largest part of the foot to pass through the smallest and tightest part of the stocking. While difficult at first, it does get easier with practice and time if the correct steps are followed.

To Apply Stockings



1 Turn leg of the stocking inside out and down as far as the heel. Except for the foot, the inside of the stocking should face outwards.



2 Use both thumbs to stretch the foot of the stocking open and pull the stocking over the foot.



3 Grip the folded edge of the stocking and pull it over the heel.



4 Once all the leg of the stocking is above the ankle, reach inside with thumbs and ease the fabric up the leg in a zig zag movement towards the knee (do not pull top seam up as this will damage the fabric).


Next, check for the following:

- For knee high stockings, there should be a two-finger width below the back of the knee and the top of the stocking.
- The heel section should fit correctly over the heel.
- The tops of the stockings should not be folded down or rolled over.
- The stocking fabric should be evenly spread along the leg to prevent fluid (from swelling or fluid retention) collecting in the loose pockets of the stockings.
- The stocking edge is positioned at the base of the toes in open toe stockings.

To Remove Stockings




1 Take hold of the top of the stocking and pull it down to the ankle.



2 Insert fingers(s) between the stocking and leg and stretch it over the ankle.



3 Slowly pull the stocking off the foot.



4 If the stocking is difficult to remove, consider a stocking aid such as an applicator, donner or glider.

See page 46.


Tips

- Apply stockings early in the morning as legs can begin to swell within 30 minutes of getting out of bed.
- Ensure leg is clean and dry.
- For those with diminished arm or hand strength, or impaired mobility there are aids to make the task of applying stockings easier. These aids are called applicators, donners, gliders and donning gloves which are all available from Independence Australia.
- Roll on stocking adhesives can be used to ensure stockings do not slide down the leg during movement.
- If not using a stocking aid, consider using a pair of kitchen style rubber gloves. Gloves can help grip the material more easily and prevent damage to the stockings and skin from sharp fingernails or rings.
- Replace stockings every 3 to 6 months, or when they start to sag or become loose.
- If new pain, discomfort, numbness or tingling is experienced in the legs or feet when wearing stockings, remove them immediately and consult your health care professional.

To access a wide range of compression stockings and stocking aids visit
store.independenceaustralia.com for quick and reliable home delivery

© Independence Australia Group 2021. No part of this information may be reproduced in whole or in part in any manner without permission.

we've got you covered!



Online Ordering
Free delivery on
orders over \$50

store.independenceaustralia.com
customerservice@iagroup.org.au
1300 788 855

your
Notes

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
BELOW KNEE									
VENOSAN 5001 STOCKING AD									
19445085	Venosan 5001 Stocking Ad	SN45011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19445105	Venosan 5001 Stocking Ad	SN45013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19445115	Venosan 5001 Stocking Ad	SN45014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19445045	Venosan 5001 Stocking Ad	SM45011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Marokko	1	1
19445065	Venosan 5001 Stocking Ad	SM45013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19445075	Venosan 5001 Stocking Ad	SM45014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19445005	Venosan 5001 Stocking Ad	SG45011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Mexico	1	1
19445025	Venosan 5001 Stocking Ad	SG45013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Mexico	1	1
19445100	Venosan 5001 Stocking Ad	SN45003	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19445040	Venosan 5001 Stocking Ad	SM45001	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Marokko	1	1
19445070	Venosan 5001 Stocking Ad	SM45004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Marokko	1	1
19445000	Venosan 5001 Stocking Ad	SG45001	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Mexico	1	1
19445020	Venosan 5001 Stocking Ad	SG45003	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Mexico	1	1
VENOSAN 5002 STOCKING AD									
19445520	Venosan 5002 Stocking Ad	SN46011	Small	23-32mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19445530	Venosan 5002 Stocking Ad	SN46013	Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19445465	Venosan 5002 Stocking Ad	SM46012	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19445470	Venosan 5002 Stocking Ad	SM46013	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19445450	Venosan 5002 Stocking Ad	SM46013	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19445475	Venosan 5002 Stocking Ad	SM46014	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19445445	Venosan 5002 Stocking Ad	SG46012	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Mexico	1	1
19445505	Venosan 5002 Stocking Ad	SN46002	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19445510	Venosan 5002 Stocking Ad	SN46003	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19445425	Venosan 5002 Stocking Ad	SM46002	Medium	34-46mmHG	Ankle 23-26cm	Open Toe	Marokko	1	1
19445435	Venosan 5002 Stocking Ad	SM46004	Extra Large	23-32mmHG	Ankle 29-31cm	Open Toe	Marokko	1	1
19445400	Venosan 5002 Stocking Ad	SG46001	Small	23-32mmHG	Ankle 29-31cm	Open Toe	Mexico	1	1
19445405	Venosan 5002 Stocking Ad	SG46002	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
19445410	Venosan 5002 Stocking Ad	SG46003	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Mexico	1	1
IMPREGNATED STOCKING									
ZIPZOC ZINC OXIDE STOCKING									
19364010	ZipZoc Zinc Oxide Stocking	66000748			80cm			1	10
KNEE HIGH									
VENOSAN 5001 STOCKING AD									
19445095	Venosan 5001 Stocking Ad	SN45012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19445055	Venosan 5001 Stocking Ad	SM45012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19445015	Venosan 5001 Stocking Ad	SG45012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Mexico	1	1
19445050	Venosan 5001 Stocking Ad	SM45.002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Marokko	1	1
19445010	Venosan 5001 Stocking Ad	SG45002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
VENOSAN 5003 STOCKING AD									
19446100	Venosan 5003 Stocking Ad	SM47002	Medium	34-46mmHG	Ankle 23-26cm	Open Toe	Marokko	1	1
19446125	Venosan 5003 Stocking Ad	SM47012	Medium	34-46mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
WAIST LENGTH									
VENOSAN 5001 STOCKING AGG									
19445225	Venosan 5001 Stocking Agg	SG45152	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
19445230	Venosan 5001 Stocking Agg	SG45172	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1



STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
KNEE HIGH								
BREVET								
19290700	Brevet Compression Stocking	8805	Medium	35-40.5cm	Open Toe	White	10	1
COMPRINET PRO ANTI EMBOLISM								
19363040	Comprinet Pro Anti Embolism Stocking	77026-00	Size 1 / Small	Calf 23-30cm	Open Toe	White	1	10
19363045	Comprinet Pro Anti Embolism Stocking	77027-00	Size 2 / Medium	Calf 30-38cm	Open Toe	White	1	10
19363050	Comprinet Pro Anti Embolism Stocking	77028-00	Size 3 / Large	Calf 38-46cm	Open Toe	White	1	10
19363055	Comprinet Pro Anti Embolism Stocking	77124-00	Size 4 / Extra Large	Calf 46-53cm	Open Toe	White	1	10
TED STOCK ANTI EMBOLISM								
19000137	Ted Stock Anti Embolism	5045	Extra Large		Close Toe	Beige	1	1
THROMBEXIN STOCK ANTI EMBOLISM								
19290010	Thrombexin Stock Anti Embolism	002-8040202	Small	Ankle 20-22cm		White	1	10
19290020	Thrombexin Stock Anti Embolism	002-8040203	Medium	Ankle 23-25cm		White	1	10
19290030	Thrombexin Stock Anti Embolism	002-8040204	Large	Ankle 26-28cm, Calf 38-48cm		White	1	10
19290035	Thrombexin Stock Anti Embolism	002-8040214	Large - Extra Wide	Ankle 29-31cm		White	1	10
19290040	Thrombexin Stock Anti Embolism	002-8040205	Extra Large			Red	1	10
LONG KNEE								
OAPL STOCK ANTI EMBOLISM								
19000715	Oapl Stock Anti Embolism	41031D	Medium	Calf 30-37cm	Open Toe	White	1	1
19000725	Oapl Stock Anti Embolism	41031F	Large	Calf 37-44cm	Open Toe	White	1	1
TED STOCK ANTI EMBOLISM								
19000126	Ted Stock Anti Embolism	4039	Large	Calf 38.2-44.5cm	Closed Toe	Beige	1	1
19000395	Ted Stock Anti Embolism	4049	Large	Calf 38.2-44.5cm	Closed Toe	Black	1	1
19000136	Ted Stock Anti Embolism	5055	Extra Large	Calf 44.6-50.8cm	Closed Toe	Black	1	1
19000105	Ted Stock Anti Embolism	7339	Small	Calf < 30.5cm	Open Toe	White	1	12
19000115	Ted Stock Anti Embolism	7480	Medium	Calf 30.5-38.1cm	Open Toe	White	1	12
19000125	Ted Stock Anti Embolism	7594	Large	Calf 38.2-44.5cm	Open Toe	White	1	12
19000135	Ted Stock Anti Embolism	7802	Extra Large	Calf 44.6-50.8cm	Open Toe	White	1	12
19000140	Ted Stock Anti Embolism	7473	XXX Large		Open Toe	White	1	12
REGULAR KNEE								
OAPL STOCK ANTI EMBOLISM								
19000700	Oapl Stock Anti Embolism	41031A	Small	Calf 25-30cm	Open Toe	White	1	1
19000710	Oapl Stock Anti Embolism	41031C	Medium	Calf 30-37cm	Open Toe	White	1	1
19000720	Oapl Stock Anti Embolism	41031E	Large	Calf 37-44cm	Open Toe	White	1	1
19000730	Oapl Stock Anti Embolism	41031G	Extra Large	Calf Over 44cm	Open Toe	White	1	1
TED STOCK ANTI EMBOLISM								
19000114	Ted Stock Anti Embolism	1087K	Medium	Calf 30.5-38.1cm	Closed Toe	Beige	1	1
19000121	Ted Stock Anti Embolism	2045	Large	Calf 38.2-44.5cm	Closed Toe	Beige	1	1
19000113	Ted Stock Anti Embolism	1097	Medium	Calf 30.5-38.1cm	Closed Toe	Black	1	1
19000122	Ted Stock Anti Embolism	2055	Large	Calf 38.2-44.5cm	Closed Toe	Black	1	1
19000132	Ted Stock Anti Embolism	3015	Extra Large	Calf 44.6-50.8cm	Closed Toe	Black	1	1
19000112	Ted Stock Anti Embolism	1077	Medium	Calf 30.5-38.1cm	Closed Toe	White	1	1
19000139	Ted Stock Anti Embolism	7470LF	XX Large		Closed Toe	White	1	12
19000102	Ted Stock Anti Embolism	1066	Small	Calf < 30.5 cm	Open Toe	Beige	1	1
19000101	Ted Stock Anti Embolism	1076	Small	Calf < 30.5 cm	Open Toe	Black	1	1
19000100	Ted Stock Anti Embolism	7071	Small	Calf < 30.5 cm	Open Toe	White	1	12
19000110	Ted Stock Anti Embolism	7115	Medium	Calf 30.5-38.1cm	Open Toe	White	1	12
19000120	Ted Stock Anti Embolism	7203	Large	Calf 38.2-44.5cm	Open Toe	White	1	12
19000130	Ted Stock Anti Embolism	7604	Extra Large	Calf 44.6-50.8cm	Open Toe	White	1	12

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
THIGH LENGTH								
COMPRINET PRO ANTI EMBOLISM								
19363070	Comprinet Pro Anti Embolism Stocking	46337-00	Size 3 / Medium	Calf 30-38cm	Open Toe	White	1	10
19363071	Comprinet Pro Anti Embolism Stocking	46338-00	Size 4 / Medium	Calf 30-38cm	Open Toe	White	1	10
19363072	Comprinet Pro Anti Embolism Stocking	46388-00	Size 5 / Large	Calf 38-46cm	Open Toe	White	1	10
19363073	Comprinet Pro Anti Embolism Stocking	46389-00	Size 6 / Large	Calf 38-46cm	Open Toe	White	1	10
19363074	Comprinet Pro Anti Embolism Stocking	77125-00	Size 7 / Extra Large	Calf 46-53cm	Open Toe	White	1	10
19363075	Comprinet Pro Anti Embolism Stocking	77128-00	Size 8 / Extra Large	Calf 46-53cm	Open Toe	White	1	10
19363060	Comprinet Pro Anti Embolism Stocking	46335-00	Small	Calf 23-30cm	Open Toe	White	1	10
TED STOCK ANTI EMBOLISM								
19000390	Ted Stock Anti Embolism	3181			Closed Toe	White	1	6
19000295	Ted Stock Anti Embolism	5707	Large	Calf 38.1-44.5cm	Closed Toe	White	1	1
19000260	Ted Stock Anti Embolism	3222LF	Small		Open Toe	White	1	6
19000255	Ted Stock Anti Embolism	3130LF	Small		Open Toe	White	1	6
19000250	Ted Stock Anti Embolism	3071LF	Small		Open Toe	White	1	6
19000270	Ted Stock Anti Embolism	3416LF	Medium		Open Toe	White	1	6
19000275	Ted Stock Anti Embolism	3549LF	Medium		Open Toe	White	1	6
19000265	Ted Stock Anti Embolism	3310LF	Medium		Open Toe	White	1	6
19000375	Ted Stock Anti Embolism	5417	Medium	Calf 30.5-38.1cm	Open Toe	White	1	1
19000290	Ted Stock Anti Embolism	3856LF	Large		Open Toe	White	1	6
19000285	Ted Stock Anti Embolism	3728LF	Large		Open Toe	White	1	6
19000280	Ted Stock Anti Embolism	3634LF	Large		Open Toe	White	1	6
THIGH LENGTH WITH BELT								
TED STOCK ANTI EMBOLISM								
19000155	Ted Stock Anti Embolism	3320	Extra Small	Calf <25.4cm	Open Toe	White	1	6
19000165	Ted Stock Anti Embolism	3364	Small	Calf 25.4-30.5cm	Open Toe	White	1	6
19000160	Ted Stock Anti Embolism	3039	Small	Calf 25.4-30.5cm	Open Toe	White	1	6
19000175	Ted Stock Anti Embolism	3449/3449P	Medium	Calf 30.5-38.1cm	Open Toe	White	1	6
19000170	Ted Stock Anti Embolism	3144	Medium	Calf 30.5-38.1cm	Open Toe	White	1	6
19000185	Ted Stock Anti Embolism	3523	Large	Calf 38.1-44.5cm	Open Toe	White	1	6
19000180	Ted Stock Anti Embolism	3221	Large	Calf 38.1-44.5cm	Open Toe	White	1	6
19000195	Ted Stock Anti Embolism	3995	Extra Large	Calf 44.5-54.6 cm	Open Toe	White	1	6
19000190	Ted Stock Anti Embolism	3922	Extra Large	Calf 44.5-54.6 cm	Open Toe	White	1	6
THROMBEXIN STOCK ANTI EMBOLISM								
19290200	Thrombexin Stock Anti Embolism	002-8062002	Small	Ankle 20-22cm		White	1	10
19290210	Thrombexin Stock Anti Embolism	002-8062003	Medium	Ankle 23-25cm		White	1	10
19290220	Thrombexin Stock Anti Embolism	002-8062004	Large	Ankle 26-28cm		White	1	5
THIGH WITH NO BELT								
THROMBEXIN STOCK ANTI EMBOLISM								
19290105	Thrombexin Stock Anti Embolism	002-8060202	Small	Ankle 20-22cm		White	1	10
19290115	Thrombexin Stock Anti Embolism	002-8060203	Medium	Ankle 23-25cm, Thigh 56-70cm		White	1	10
19290120	Thrombexin Stock Anti Embolism	0028060213	Medium - Extra Wide	Ankle 23-25cm, Thigh 64-80cm		White	1	10
19290125	Thrombexin Stock Anti Embolism	002-8060204	Large	Ankle 26-28cm, Thigh 64-80cm		White	1	1
19290130	Thrombexin Stock Anti Embolism	002-8060214	Large - Extra Wide	Ankle 26-28cm, Thigh 72-90cm		White	1	10
19290135	Thrombexin Stock Anti Embolism	002-8060205	Extra Large	Ankle 29-31cm		White	1	10

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
ABDOMINAL BINDER									
ABDOMINAL BINDER									
19310110	Abdominal Binder	SM829	Small		Height 22.5cm		White	1	1
19310100	Abdominal Binder	SM830	Regular		Height 15cm		White	1	1
19310120	Abdominal Binder	SM832	Regular		Height 22.5cm		White	1	1
19310140	Abdominal Binder	SM834	Regular		Height 30cm		White	1	1
19310105	Abdominal Binder	SM831	Large		Height 15cm		White	1	1
19310125	Abdominal Binder	SM833	Large		Height 22.5cm		White	1	1
19310145	Abdominal Binder	SM835	Large		Height 30cm		White	1	1
19310130	Abdominal Binder	SM837	Extra Large		Height 22.5cm		White	1	1
19310146	Abdominal Binder	SM836	Extra Large		Height 30cm		White	1	1
ARMSLEEVE									
JOBST									
19360045	Jobst Bella Armsleeve	72737-05	Size 5	20-30mmHG	Palm 19-21cm Wrist 18-21cm			1	1
19363085	Jobst RTW Sleeve	75277-00	Medium	20-30mmHG	Mid Upper 27-45cm Wrist 16-22cm			1	1
19363090	Jobst RTW Sleeve	75278-00	Large	20-30mmHG	Mid Upper 33-51cm Wrist 18-24cm			1	1
ARMSLEEVE - CUSTOM MADE PRODUCT									
CAROSKIN ARM SLEEVE									
19000005	Caroskin Arm Sleeve Class 1	4012001	Custom Made	18-21mmHG				1	1
19000015	Caroskin Arm Sleeve Class 1	4010001	Custom Made	18-21mmHG				1	1
19000030	Caroskin Arm Sleeve Class 2	4012001	Custom Made	18-21mmHG				1	1
19000040	Caroskin Arm Sleeve Class 3	4012001	Custom Made	34-46mmHG				1	1
COMPRESSION GLOVES									
VENOSAN COMPRESSION GLOVES									
19448520	Venosan Compression Gloves 3/4 Finger	VE94521	Small Plus	23-32mmHG	Mild 16-19cm		Beige	1	1
19448525	Venosan Compression Gloves 3/4 Finger	VE94522	Medium	23-32mmHG	Mild 16-19cm		Beige	1	1
COMPRESSION GLOVES - CUSTOM MADE PRODUCT									
CAROSKIN									
19000020	Caroskin Glove With No Fingers Class 1	4016001	Custom Made	18-21mmHG				1	1
19000025	Caroskin Gauntlet No Fingers Class 1	4017001	Custom Made	18-21mmHG				1	1
GRIP TOP STOCKING - CUSTOM MADE PRODUCT									
CAROSKIN GRIP TOP									
19000010	Caroskin Grip Top Class 1	5018003	Custom Made		Band 3cm			1	1
19000035	Caroskin Grip Top Class 1	5018003	Custom Made		Band 3cm			1	1
19000045	Caroskin Grip Top Class 1	5018003	Custom Made		Band 3cm	Open Toe		1	1
OEDEMA CONTROL - COMPRESSION GARMENTS									
FARROW WRAP LITE									
19000220	Farrow Wrap Lite Leg Piece	76654-01	Small	20-30mmHG				1	1
19000225	Farrow Wrap Lite Foot Piece	76652-04	Extra Large	20-30mmHG				1	1

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
SLIP RESISTANT SOCKS									
SOCKS SLIP RESISTANT									
19008455	Socks Slip Resistant	NSSBLCTN20/FPSH-SRS-LGBU	Large				Burgundy	20	20
19008460	Socks Slip Resistant	ST NSSRLCTN20	Large		Size 10-13		Red	1	1
19008450	Socks Slip Resistant Economy	FPSH-ESRS-XLRDCTN48/ FPSH-ESRS-OSRDPR	Extra Large				Red	1	1
SUPPORT SOCKS									
JOBST COMPRESSION									
19360000	Jobst Compression Socks	75253-01	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360100	Jobst Compression Socks	75257-01	Small	30-40mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360200	Jobst Compression Socks	75265-01	Small	30-40mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360350	Jobst Compression Socks	75261-00	Small	30-40mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360005	Jobst Compression Socks	75254-01	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360105	Jobst Compression Socks	75258-01	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360205	Jobst Compression Socks	75266-01	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360285	Jobst Compression Socks	75270-00	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360010	Jobst Compression Socks	75255-01	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360110	Jobst Compression Socks	75259-01	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360210	Jobst Compression Socks	75267-01	Large	30-40mmHG	Ankle 23-29cm	Closed Toe	Black	1	1
19360290	Jobst Compression Socks	75271-00	Large	30-40mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360015	Jobst Compression Socks	75256-01	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360115	Jobst Compression Socks	75260-01	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360295	Jobst Compression Socks	75272-00	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360215	Jobst Compression Socks	75268-01	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
VENOSAN MICROFIBRE									
19444290	Venosan Microfibre Socks	MC79205	Small	15-20mmHG	Ankle 19-21cm	Closed Toe	Beige	1	1
19444295	Venosan Microfibre Socks	MC79206	Medium	15-20mmHG	Ankle 21-24cm	Closed Toe	Beige	1	1
19444300	Venosan Microfibre Socks	MC79207	Large	15-20mmHG	Ankle 24-26cm	Closed Toe	Beige	1	1
19444305	Venosan Microfibre Socks	MC79208	Extra Large	15-20mmHG	Ankle 26-32cm	Closed Toe	Beige	1	1
19444310	Venosan Microfibre Socks	MC79221	Small	15-20mmHG	Ankle 19-21cm	Closed Toe	Black	1	1
19444315	Venosan Microfibre Socks	MC79222	Medium	15-20mmHG	Ankle 21-24cm	Closed Toe	Black	1	1
19444255	Venosan Microfibre Socks	MC79003	Medium	15-20mmHG	Ankle 21-24cm	Closed Toe	Black	1	1
19444275	Venosan Microfibre Socks	MC89003	Medium	20-30mmHG	Ankle 21-24cm	Closed Toe	Black	1	1
19444320	Venosan Microfibre Socks	MC79223	Large	15-20mmHG	Ankle 24-26cm	Closed Toe	Black	1	1
19444260	Venosan Microfibre Socks	MC79004	Large	15-20mmHG	Ankle 22-26cm	Closed Toe	Black	1	1
19444280	Venosan Microfibre Socks	MC89004	Large	20-30mmHG	Ankle 22-26cm	Closed Toe	Black	1	1
19444285	Venosan Microfibre Socks	MC89005	Extra Large	20-30mmHG	Ankle 26-32cm	Closed Toe	Black	1	1
VENOSAN SILVERLINE									
19444150	Venosan Silverline Socks	VS49511	Small	20-30mmHG	Ankle 19-21cm	Closed Toe	Beige	1	1
19444155	Venosan Silverline Socks	VS49512	Medium	20-30mmHG	Ankle 21-24cm	Closed Toe	Beige	1	1
19444160	Venosan Silverline Socks	VS49513	Large	30-40mmHG	Ankle 24-26cm	Closed Toe	Beige	1	1
19444165	Venosan Silverline Socks	VS49514	Extra Large	30-40mmHG	Ankle 26-32cm	Closed Toe	Beige	1	1
19444170	Venosan Silverline Socks	VS49515	XX Large	20-30mmHG	Ankle 32-36cm	Closed Toe	Beige	1	1
19444100	Venosan Silverline Socks	VS49405	Small	20-30mmHG	Ankle 19-21cm	Closed Toe	Black	1	1
19444130	Venosan Silverline Socks	VS49506	Medium	20-30mmHG	Ankle 21-24cm	Closed Toe	Black	1	1
19444105	Venosan Silverline Socks	VS49406	Medium	20-30mmHG	Ankle 21-24cm	Closed Toe	Black	1	1
19444135	Venosan Silverline Socks	VS49507	Large	30-40mmHG	Ankle 24-26cm	Closed Toe	Black	1	1
19444110	Venosan Silverline Socks	VS49407	Large	30-40mmHG	Ankle 24-26cm	Closed Toe	Black	1	1
19444140	Venosan Silverline Socks	VS49508	Extra Large	30-40mmHG	Ankle 26-32cm	Closed Toe	Black	1	1
19444115	Venosan Silverline Socks	VS49408	Extra Large	30-40mmHG	Ankle 26-32cm	Closed Toe	Black	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Stock Code	Product Description	Supp. Part	Size	Compression	Measurement	Style	Colour	Packet Qty	Carton Qty
19444120	Venosan Silverline Socks	VS49409	XX Large	20-30mmHG	Ankle 32-36cm	Closed Toe	Black	1	1
VENOSAN SPORTSLINE									
19446700	Venosan Sportslinesocks	SP69202	Small	20-30mmHG	Ankle 19-23cm	Closed Toe	White	1	1
19446705	Venosan Sportslinesocks	SP69203	Medium	15-20mmHG	Ankle 23-26cm	Closed Toe	White	1	1
19446725	Venosan Sportslinesocks	SP79203	Medium	20-30mmHG	Ankle 23-26cm	Closed Toe	White	1	1
19446710	Venosan Sportslinesocks	SP69204	Large	15-20mmHG	Ankle 26-29cm	Closed Toe	White	1	1
19446730	Venosan Sportslinesocks	SP79204	Large	20-30mmHG	Ankle 26-29cm	Closed Toe	White	1	1
19446715	Venosan Sportslinesocks	SP69205	Extra Large	15-20mmHG	Ankle 29-31cm	Closed Toe	White	1	1
TRAVEL SOCKS									
JOBST ACTIVE WEAR SOCKS									
19362950	Jobst Active Wear Socks	75156-03	Medium	20-30mmHG		Closed Toe	Black	1	1
JOBST TRAVEL									
19363660	Jobst Travel Socks	78844-07	Size 1	15-20mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19363645	Jobst Travel Socks	78845-12	Size 2	15-20mmHG	Ankle 20-23cm	Closed Toe	Beige	1	1
19363665	Jobst Travel Socks	78845-07	Size 2	15-20mmHG	Ankle 20-23cm	Closed Toe	Black	1	1
19363650	Jobst Travel Socks	78846-12	Size 3	15-20mmHG	Ankle 22-25cm	Closed Toe	Beige	1	1
19363670	Jobst Travel Socks	78846-07	Size 3	15-20mmHG	Ankle 22-25cm	Closed Toe	Black	1	1
19363655	Jobst Travel Socks	78847-12	Size 4	15-20mmHG	Ankle 24-27cm	Closed Toe	Beige	1	1
19363675	Jobst Travel Socks	78847-07	Size 4	15-20mmHG	Ankle 24-27cm	Closed Toe	Black	1	1
19363680	Jobst Travel Socks	79272-07	Size 5	15-20mmHG	Ankle 26-29cm	Closed Toe	Black	1	1

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
BELOW KNEE									
OAPL THERMAFIRM									
19000685	Oapl Thermafirm Stocking	68202	Large	15-20mmHG	Ankle 23-27cm	Closed Toe	Sand	1	1
19000605	Oapl Thermafirm Stocking	75104/75113	Medium	20-30mmHG		Closed Toe	White	1	1
19000670	Oapl Thermafirm Stocking		XX Large	30-40mmHG		Open Toe	Beige	1	1
VENOSAN 4001									
19442780	Venosan 4001 Stocking Ad	SN41011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19442785	Venosan 4001 Stocking Ad	SN41012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19442790	Venosan 4001 Stocking Ad	SN41013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19442795	Venosan 4001 Stocking Ad	SN41014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19442660	Venosan 4001 Stocking Ad	SM41011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Marokko	1	1
19442665	Venosan 4001 Stocking Ad	SM41012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19442670	Venosan 4001 Stocking Ad	SM41013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19442675	Venosan 4001 Stocking Ad	SM41014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19442605	Venosan 4001 Stocking Ad	SG41011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Mexico	1	1
19442615	Venosan 4001 Stocking Ad	SG41012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Mexico	1	1
19442625	Venosan 4001 Stocking Ad	SG41013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Mexico	1	1
19442760	Venosan 4001 Stocking Ad	SN41001	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Black	1	1
19442765	Venosan 4001 Stocking Ad	SN41002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19442770	Venosan 4001 Stocking Ad	SN41003	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19442775	Venosan 4001 Stocking Ad	SN41004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Black	1	1
19442640	Venosan 4001 Stocking Ad	SM41001	Small	18-21mmHG	Ankle19-23cm	Open Toe	Marokko	1	1
19442645	Venosan 4001 Stocking Ad	SM41002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Marokko	1	1
19442650	Venosan 4001 Stocking Ad	SM41003	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Marokko	1	1
19442655	Venosan 4001 Stocking Ad	SM41004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Marokko	1	1
19442600	Venosan 4001 Stocking Ad	SG41001	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Mexico	1	1
19442610	Venosan 4001 Stocking Ad	SG41002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
19442620	Venosan 4001 Stocking Ad	SG41003	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Mexico	1	1
19442630	Venosan 4001 Stocking Ad	SG41004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Mexico	1	1
VENOSAN 4002									
19442910	Venosan 4002 Stocking Ad	SN42011	Small	25-32mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19442960	Venosan 4002 Stocking Ad	SN42012	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19443010	Venosan 4002 Stocking Ad	SN42013	Large	25-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19443060	Venosan 4002 Stocking Ad	SN42014	Extra Large	25-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19442880	Venosan 4002 Stocking Ad	SM42011	Small	25-32mmHG	Ankle 19-23cm	Closed Toe	Marokko	1	1
19442930	Venosan 4002 Stocking Ad	SM42012	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19442980	Venosan 4002 Stocking Ad	SM42013	Large	25-32mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19443030	Venosan 4002 Stocking Ad	SM42014	Extra Large	25-32mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19442870	Venosan 4002 Stocking Ad	SG42011	Small	25-32mmHG	Ankle 19-23cm	Closed Toe	Mexico	1	1
19442920	Venosan 4002 Stocking Ad	SG42012	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Mexico	1	1
19442970	Venosan 4002 Stocking Ad	SG42013	Large	25-32mmHG	Ankle 26-29cm	Closed Toe	Mexico	1	1
19443020	Venosan 4002 Stocking Ad	SG42014	Extra Large	25-32mmHG	Ankle 29-31cm	Closed Toe	Mexico	1	1
19442955	Venosan 4002 Stocking Ad	SN42002	Medium	25-32mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19443005	Venosan 4002 Stocking Ad	SN42003	Large	25-32mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19443055	Venosan 4002 Stocking Ad	SN42004	Extra Large	25-32mmHG	Ankle 29-31cm	Open Toe	Black	1	1
19442875	Venosan 4002 Stocking Ad	SM42001	Small	25-32mmHG	Ankle 19-23cm	Open Toe	Marokko	1	1
19442925	Venosan 4002 Stocking Ad	SM42002	Medium	25-32mmHG	Ankle 23-26cm	Open Toe	Marokko	1	1
19442975	Venosan 4002 Stocking Ad	SM42.003	Large	25-32mmHG	Ankle 26-29cm	Open Toe	Marokko	1	1
19443025	Venosan 4002 Stocking Ad	SM42004	Extra Large	25-32mmHG	Ankle 29-31cm	Open Toe	Marokko	1	1
19442865	Venosan 4002 Stocking Ad	SG42001	Small	25-32mmHG	Ankle 19-23cm	Open Toe	Mexico	1	1
19442915	Venosan 4002 Stocking Ad	SG42002	Medium	25-32mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
19442965	Venosan 4002 Stocking Ad	SG42003	Large	25-32mmHG	Ankle 26-29cm	Open Toe	Mexico	1	1
19443015	Venosan 4002 Stocking Ad	SG42004	Extra Large	25-32mmHG	Ankle 29-31cm	Open Toe	Mexico	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
VENOSAN LEGLINE									
19446265	Venosan Legline 20 Ad	VL406BL	Large	30mmHG	Ankle 26-27cm	Closed Toe	Black	1	1
19446015	Venosan Legline 20 Ad	VL304BL	Small	20mmHG	Ankle 19-23cm	Closed Toe Sheer	Black	1	1
19446226	Venosan Legline 30 Ad	VL405BL	Medium	30mmHG	Ankle 24-25cm	Closed Toe Sheer	Black	1	1
19446227	Venosan Legline 30 Ad	VL406BL	Large	30mmHG	Ankle 26-27cm	Closed Toe Sheer	Black	1	1
19446220	Venosan Legline 30 Ad	VL403NU	Extra Small	30mmHG	Ankle 20-21cm	Closed Toe Sheer	Nude	1	1
19446045	Venosan Legline 20 Ad	VL305BL	Medium	20mmHG	Ankle 24-25cm	Closed Toe Sheer	Nude	1	1
19446070	Venosan Legline 20 Ad	VL305NU	Medium	20mmHG	Ankle 24-25cm	Closed Toe Sheer	Nude	1	1
19446221	Venosan Legline 30 Ad	VL404NU	Medium	30mmHG	Ankle 22-23cm	Closed Toe Sheer	Nude	1	1
19446075	Venosan Legline 20 Ad	VL306NU	Large	20mmHG	Ankle 26-27cm	Closed Toe Sheer	Nude	1	1
19446223	Venosan Legline 30 Ad	VL406NU	Large	30mmHG	Ankle 26-27cm	Closed Toe Sheer	Nude	1	1
19446020	Venosan Legline 20 Ad	VL304SA	Small	20mmHG	Ankle 22-23cm	Closed Toe Sheer	Sahara	1	1
19446024	Venosan Legline 20 Ad	VL305SA	Medium	20mmHG	Ankle 24-25cm	Closed Toe Sheer	Sahara	1	1
19446230	Venosan Legline 30 Ad	VL405SA	Medium	30mmHG	Ankle 24-25cm	Closed Toe Sheer	Sahara	1	1
19446231	Venosan Legline 30 Ad	VL406SA	Large	30mmHG	Ankle 26-27cm	Closed Toe Sheer	Sahara	1	1
19446025	Venosan Legline 20 Ad	VL306SA	Large	20mmHG	Ankle 26-27cm	Closed Toe Sheer	Sahara	1	1
VENOSAN 6001									
19446360	Venosan 6001 Stocking Ad	CG61031	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Beige	1	1
19446320	Venosan 6001 Stocking Ad	CG61011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Beige	1	1
19446325	Venosan 6001 Stocking Ad	CG61012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Beige	1	1
19446369	Venosan 6001 Stocking Ad	CG61032	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Beige	1	1
19446370	Venosan 6001 Stocking Ad	CG61033	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19446375	Venosan 6001 Stocking Ad	CG61034	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
19446335	Venosan 6001 Stocking Ad	CG61014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
19446400	Venosan 6001 Stocking Ad	CN61011	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19446405	Venosan 6001 Stocking Ad	CN61012	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19446445	Venosan 6001 Stocking Ad	CN61032	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19446455	Venosan 6001 Stocking Ad	CN61034	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19446415	Venosan 6001 Stocking Ad	CN61014	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19446300	Venosan 6001 Stocking Ad	CG61001	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19446340	Venosan 6001 Stocking Ad	CG61021	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19446345	Venosan 6001 Stocking Ad	CG61022	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446305	Venosan 6001 Stocking Ad	CG61002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446315	Venosan 6001 Stocking Ad	CG61004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Beige	1	1
19446355	Venosan 6001 Stocking Ad	CG61024	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Beige	1	1
19446420	Venosan 6001 Stocking Ad	CN61021	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Black	1	1
19446385	Venosan 6001 Stocking Ad	CN61002	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19446395	Venosan 6001 Stocking Ad	CN61004	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Black	1	1
19446435	Venosan 6001 Stocking Ad	CN61024	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Black	1	1
VENOSAN 6002									
19446600	Venosan 6002 Stocking Ad	CN62031	Small	23-32mmHG	Ankle 19-23cm	Closed Toe	Black	1	1
19446585	Venosan 6002 Stocking Ad	CN62012	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19446605	Venosan 6002 Stocking Ad	CN62032	Medium	23-32mmHG	Ankle 23-27cm	Closed Toe	Black	1	1
19446595	Venosan 6002 Stocking Ad	CN62014	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19446615	Venosan 6002 Stocking Ad	CN62034	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19446500	Venosan 6002 Stocking Ad	CG62011	Small	23-32mmHG	Ankle 19-23cm	Closed Toe	Beige	1	1
19446520	Venosan 6002 Stocking Ad	CG62031	Small	23-32mmHG	Ankle 19-23cm	Closed Toe	Beige	1	1
19446505	Venosan 6002 Stocking Ad	CG62012	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Beige	1	1
19446525	Venosan 6002 Stocking Ad	CG62032	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Beige	1	1
19446535	Venosan 6002 Stocking Ad	CG62034	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
19446515	Venosan 6002 Stocking Ad	CG62014	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
19446540	Venosan 6002 Stocking Ad	CN62001	Small	23-32mmHG	Ankle 19-23cm	Open Toe	Black	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
19446545	Venosan 6002 Stocking Ad	CN62002	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19446555	Venosan 6002 Stocking Ad	CN62004	Extra Large	23-32mmHG	Ankle 29-31cm	Open Toe	Black	1	1
19446575	Venosan 6002 Stocking Ad	CN62024	Extra Large	23-32mmHG	Ankle 29-31cm	Open Toe	Black	1	1
19446460	Venosan 6002 Stocking Ad	CG62001	Small	23-32mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19446465	Venosan 6002 Stocking Ad	CG62002	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446485	Venosan 6002 Stocking Ad	CG62022	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446475	Venosan 6002 Stocking Ad	CG62004	Extra Large	23-32mmHG	Ankle 29-31cm	Open Toe	Beige	1	1
VENOSAN 6003									
19446625	Venosan 6003 Stocking Ad	CG63002	Medium	34-46mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446655	Venosan 6003 Stocking Ad	CG63024	Extra Large	34-46mmHG	Ankle 29-31cm	Open Toe	Beige	1	1
VENOSAN 7001									
19447005	Venosan 7001 Stocking	SG71021	Small	18-21mmHG	Ankle 20-23cm	Open Toe	Beige	1	1
VENOSAN LEGLINE									
19447220	Venosan Legline	VL336BL	Large	20mmHG	Ankle 26-27cm	Closed Toe	Black	1	1
19447245	Venosan Legline	VL336SA	Large	20mmHG	Ankle 26-27cm	Closed Toe	Sahara	1	1
19447250	Venosan Legline	VL337SA	Extra Large	20mmHG	Ankle 28-29cm	Closed Toe	Sahara	1	1
CALF LENGTH									
JOBST RELIEF HOSIERY									
19362011	Jobst Relief Hosiery	114698A	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Beige	1	1
19361485	Jobst Relief Hosiery	75344-12(114817)	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe (Pair)	Black	1	1
19362020	Jobst Relief Hosiery	114624A/77976-16	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Beige	1	1
19362070	Jobst Relief Hosiery	114634A	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Beige	1	1
19361505	Jobst Relief Hosiery	75343-10(114804)	Large	15-20mmHG	Ankle 25-29cm	Open Toe (Pair)	Beige	1	1
19362115	Jobst Relief Hosiery	114696A	Large	20-30mmHG	Ankle 25-29cm	Open Toe	Beige	1	1
19362125	Jobst Relief Hosiery	114629A	Extra Large	20-30mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
19362165	Jobst Relief Hosiery	114697A	Extra Large	30-40mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
19362175	Jobst Relief Hosiery	77978-15	Extra Large	30-40mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
SIGVARIS									
19001015	Sigvaris Stock Traditional	10118	Large	23-32mmHG	Calf A-D	Open Toe	Beige	1	1
19001025	Sigvaris Stock Traditional	10109	Large	23-32mmHG	Calf A-D	Open Toe	Beige	1	1
CHAPS									
JOBST RELIEF									
19362455	Jobst Relief Hosiery	114673A	Medium	20-30mmHG	Ankle 21-25cm	Open Toe	Beige	1	1
19361635	Jobst Relief Hosiery	114790A	Large	30-40mmHG		Open Toe (Pair)	Beige	1	1
19363105	Jobst Relief Chaps	114674A	Large Left	20-30mmHG	Hip 81-137cm		Beige	1	1
KNEE HIGH									
JOBST CASUAL VASCULAR									
19365280	Jobst Casual Vascular Hosiery	75485-03	Medium	20-30mmHG	Ankle 21cm-25cm	Closed Toe	Khaki	1	1
19362960	Jobst Casual Vascular Hosiery	75478-03	Large	15-20mmHG	Ankle 25cm-29cm	Closed Toe	Khaki	1	1
19362970	Jobst Casual Vascular Hosiery	75479-03	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Khaki	1	1
19362955	Jobst Casual Vascular Hosiery	75478-02	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19365270	Jobst Casual Vascular Hosiery	75486-02	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19365275	Jobst Casual Vascular Hosiery	75487-02	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19362965	Jobst Casual Vascular Hosiery	75479-02	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
JOBST ACTIVE VASCULAR									
19365115	Jobst Active Vascular Hosiery	75148-03	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19365165	Jobst Active Vascular Hosiery	75157-03	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19365205	Jobst Active Vascular Hosiery	75413-03	Large	30-40mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19365175	Jobst Active Vascular Hosiery	75158-03	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19365135	Jobst Active Vascular Hosiery	75150-03	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19365160	Jobst Active Vascular Hosiery	75157-02	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	White	1	1
19365210	Jobst Active Vascular Hosiery	75414-02	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	White	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
JOBST HOSIERY ULTRASHEER									
19360570	Jobst Hosiery Ultrasheer	75167-30	Large	15-20mmHG	Ankle 25-29cm	Open Toe	Beige	1	1
19360400	Jobst Hosiery Ultrasheer	75165-04	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360800	Jobst Hosiery Ultrasheer	75201-04	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360865	Jobst Hosiery Ultrasheer	73634-01	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360405	Jobst Hosiery Ultrasheer	75166-04	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360805	Jobst Hosiery Ultrasheer	75202-04	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19361205	Jobst Hosiery Ultrasheer	75242-04	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360410	Jobst Hosiery Ultrasheer	75167-04	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360810	Jobst Hosiery Ultrasheer	75203-04	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360415	Jobst Hosiery Ultrasheer	75168-04	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360815	Jobst Hosiery Ultrasheer	75204-04	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19361215	Jobst Hosiery Ultrasheer	75244-04	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360440	Jobst Hosiery Ultrasheer	75165-25	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19360820	Jobst Hosiery Ultrasheer	75201-19	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19360445	Jobst Hosiery Ultrasheer	75166-25	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360825	Jobst Hosiery Ultrasheer	75202-19	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19361225	Jobst Hosiery Ultrasheer	75242-21	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360450	Jobst Hosiery Ultrasheer	75167-25	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19360830	Jobst Hosiery Ultrasheer	75203-19	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19361230	Jobst Hosiery Ultrasheer	75243-21	Large	30-40mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19360455	Jobst Hosiery Ultrasheer	75168-25	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
19360835	Jobst Hosiery Ultrasheer	75204-19	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
19361235	Jobst Hosiery Ultrasheer	75244-21	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
JOBST RELIEF									
19000050	Jobst Relief KH M CT	78051-15	Medium	20-30mmHg	Ankle 21-25cm	Closed Toe	Beige	1	1
19000055	Jobst Relief KH LG CT	78052-15	Large	20-30mmHg	Ankle 25-29cm	Closed Toe	Beige	1	1
19000070	Jobst Relief KH LG OT	78052-04	Large	20-30mmHg	Ankle 25-29cm	Open Toe	Beige	1	1
19000065	Jobst Relief KH XL CT	78053-15	Extra Large	20-30mmHg	Ankle 29-33cm	Closed Toe	Beige	1	1
19000080	Jobst Relief KH XL OT	78053-04	Extra Large	20-30mmHg	Ankle 29-33cm	Open Toe	Beige	1	1
19000060	Jobst Relief KH FC LG CT	77975-16	Full Calf	20-30mmHg		Closed Toe	Beige	1	1
19000075	Jobst Relief KH FC LG OT	77975-15	Full Calf	20-30mmHg		Open Toe	Beige	1	1
19000085	Jobst Relief KH FC XL OT	77976-15	Full Calf	20-30mmHg		Open Toe	Beige	1	1
JOBST RELIEF HOSIERY									
19361430	Jobst Relief Hosiery	75339-11(114806)	Small	15-20mmHG	Ankle 18-21cm	Closed Toe (Pair)	Beige	1	1
19362000	Jobst Relief Hosiery	114620A	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Beige	1	1
19361435	Jobst Relief Hosiery	75340-11(114807)	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe (Pair)	Beige	1	1
19362005	Jobst Relief Hosiery	114621A	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Beige	1	1
19362055	Jobst Relief Hosiery	114631A	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Beige	1	1
19361440	Jobst Relief Hosiery	75341-11(114808)	Large	15-20mmHG	Ankle 25-29cm	Closed Toe (Pair)	Beige	1	1
19362010	Jobst Relief Hosiery	114622A	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Beige	1	1
19362015	Jobst Relief Hosiery	114623A	Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Beige	1	1
19362060	Jobst Relief Hosiery	114632A	Large	30-40mmHG	Ankle 25-29cm	Closed Toe	Beige	1	1
19362065	Jobst Relief Hosiery	114633A	Extra Large	30-40mmHG	Ankle 29-33cm	Closed Toe	Beige	1	1
19361460	Jobst Relief Hosiery	75339-12(114812)	Small	15-20mmHG	Ankle 18-21cm	Closed Toe (Pair)	Black	1	1
19361465	Jobst Relief Hosiery	75340-12 (114813)	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe (Pair)	Black	1	1
19361470	Jobst Relief Hosiery	75341-12(114814)	Large	15-20mmHG	Ankle 25-29cm	Closed Toe (Pair)	Black	1	1
19362100	Jobst Relief Hosiery	78050-04	Small	20-30mmHG	Ankle 18-21cm	Open Toe	Beige	1	1
19362150	Jobst Relief Hosiery	114635A	Small	30-40mmHG	Ankle 18-21cm	Open Toe	Beige	1	1
19361495	Jobst Relief Hosiery	75340-10(114801)	Medium	15-20mmHG	Ankle 21-25cm	Open Toe (Pair)	Beige	1	1
19362105	Jobst Relief Hosiery	78051-04	Medium	20-30mmHG	Ankle 21-25cm	Open Toe	Beige	1	1
19362155	Jobst Relief Hosiery	78055-04	Medium	30-40mmHG	Ankle 21-25cm	Open Toe	Beige	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
19362110	Jobst Relief Hosiery	114627A	Large	20-30mmHG	Ankle 25-29cm	Open Toe	Beige	1	1
19362160	Jobst Relief Hosiery	114637A	Large	30-40mmHG	Ankle 25-29cm	Open Toe	Beige	1	1
19361510	Jobst Relief Hosiery	75342-10(114803)	Extra Large	15-20mmHG	Ankle 29-33cm	Open Toe (Pair)	Beige	1	1
19362120	Jobst Relief Hosiery	114628A	Extra Large	20-30mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
19362170	Jobst Relief Hosiery	114638A	Extra Large	30-40mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
HEALTH SUPPORT SYSTEM									
19008000	Health Support System Hosiery	HSML 30-35 A R Beige	Size A	Class 2	Ankle 18-20cm	Closed Toe	Beige	1	1
19008010	Health Support System Hosiery	HSML 30-35 B R Beige	Size B	Class 2	Ankle 20-22.5cm	Closed Toe	Beige	1	1
19008015	Health Support System Hosiery	HSML 30-35 B S Beige	Size B	Class 2	Ankle 20-22.5cm	Closed Toe	Beige	1	1
19008020	Health Support System Hosiery	HSML 30-35 C R Beige	Size C	Class 2	Ankle 22.5-25.5cm	Closed Toe	Beige	1	1
19008030	Health Support System Hosiery	HSML 30-35 D R Beige	Size D	Class 2	Ankle 25.5-28cm	Closed Toe	Beige	1	1
19008035	Health Support System Hosiery	HSML 30-35 D S Beige	Size D	Class 2	Ankle 25.5-28cm	Closed Toe	Beige	1	1
19008040	Health Support System Hosiery	HSML 30-35 E R Beige	Size E	Class 2	Ankle 28-30.5cm	Closed Toe	Beige	1	1
19008050	Health Support System Hosiery	HSML 30-35 F R Beige	Size F	Class 2	Ankle 30.5-33cm	Closed Toe	Beige	1	1
19008160	Health Support System Hosiery	HSML 35-40 E R Beige	Size E	Class 3	Ankle 28-30.5cm	Closed Toe	Beige	1	1
19008165	Health Support System Hosiery	HSML 35-40 E S Beige	Size E	Class 3	Ankle 28-30.5cm	Closed Toe	Beige	1	1
SIGVARIS STOCK									
19000000	Sigvaris Stock Traditional 503 Class 2	10120	Extra Large Plus	23-32mmHG			Beige	1	1
19001006	Sigvaris Stock Cotton	50370	Medium	23-32mmHG	Calf A-D	Closed Toe	Black	1	1
19000970	Sigvaris Stock Cotton	63944	Medium	18-21mmHG	Calf A-D	Closed Toe	Nature	1	1
19001150	Sigvaris Stock Cotton	63945	Medium	18-21mmHG	Calf A-D	Closed Toe	Nature	1	1
19000885	Sigvaris Stock Cotton	54036	Small Plus	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19000985	Sigvaris Stock Cotton	63941	Small	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19001010	Sigvaris Stock Cotton	54037	Small	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19000875	Sigvaris Stock Cotton	50304	Small	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000940	Sigvaris Stock Cotton	50305	Small	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000935	Sigvaris Stock Cotton	54039	Medium	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19000880	Sigvaris Stock Cotton	50306	Medium	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000910	Sigvaris Stock Cotton	50307	Medium	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000930	Sigvaris Stock Cotton	50315	Medium Plus	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19001045	Sigvaris Stock Cotton	50362	Medium Plus	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19001095	Sigvaris Stock Cotton	54041	Large	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19000955	Sigvaris Stock Cotton	54040	Large	18-21mmHG	Calf A-D	Open Toe	Nature	1	1
19000900	Sigvaris Stock Cotton	50308	Large	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000915	Sigvaris Stock Cotton	50309	Large	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000965	Sigvaris Stock Cotton	50316	Large Plus	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
19000980	Sigvaris Stock Cotton	50317	Large Plus	23-32mmHG	Calf A-D	Open Toe	Nature	1	1
VENOSAN LEGLINE									
19446065	Venosan Legline 20 Ad	VL304NU	Small	20mmHG	Ankle 22-23cm	Closed Toe Sheer	Nude	1	1
VENOSAN 6001									
19446310	Venosan 6001 Stocking Ad	CG61003	Large	18-21mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19446330	Venosan 6001 Stocking Ad	CG61013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19446350	Venosan 6001 Stocking Ad	CG61023	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Beige	1	1
19446410	Venosan 6001 Stocking Ad	CN61013	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19446430	Venosan 6001 Stocking Ad	CN61023	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Black	1	1
VENOSAN 6002									
19446510	Venosan 6002 Stocking Ad	CG62013	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19446530	Venosan 6002 Stocking Ad	CG62033	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19446590	Venosan 6002 Stocking Ad	CN62013	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19446610	Venosan 6002 Stocking Ad	CN62033	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19446470	Venosan 6002 Stocking Ad	CG62003	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Beige	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
19446490	Venosan 6002 Stocking Ad	CG62023	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Beige	1	1
19446550	Venosan 6002 Stocking Ad	CN62003	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19446570	Venosan 6002 Stocking Ad	CN62023	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Black	1	1
VENOSAN 6003									
19446630	Venosan 6003 Stocking Ad	CG63003	Large	34-46mmHG	Ankle 26-29cm	Open Toe	Beige	1	1
19446650	Venosan 6003 Stocking Ad	CG63023	Large	34-46mmHG	Ankle 26-29cm	Open Toe	Beige	1	1
VENOSAN STOCKING									
19446695	Venosan Stocking Ad	CG63034	Extra Large	34-46mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
LINERS									
HEALTH SUPPORT SYSTEM									
19008260	Health Support System Liner	LUCKi W-inner Reg B	Size B		Ankle 20-22.5cm		White	1	1
MATERNITY									
JOBST HOSIERY									
19361360	Jobst Hosiery Ultrasheer	75180-01	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360700	Jobst Hosiery Ultrasheer	75177-04	Small	15-20mmHG	Ankle 18-22cm	Closed Toe	Natural	1	1
19360760	Jobst Hosiery Ultrasheer	75178-04	Medium	15-20mmHG	Ankle 21-25cm		Natural	1	1
19360755	Jobst Hosiery Ultrasheer	75179-04	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19360715	Jobst Hosiery Ultrasheer	75180-04	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
VENOSAN 4001									
19447550	Venosan 4001 Pantyhose	SG41304	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Mexico	1	1
VENOSAN 4002									
19443347	Venosan 4002 Pantyhose	SG42301	Small	23-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19443350	Venosan 4002 Pantyhose	SN42303	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
PANTYHOSE									
VENOSAN LEGLINE									
19447165	Venosan Legline 30 Panythose	VL435BL	Medium	30mmHG	Ankle 24-25cm	Closed Toe Sheer	Black	1	1
19447100	Venosan Legline 20 Panythose	VL337BL	Extra Large	20mmHG	Ankle 29-31cm	Closed Toe Sheer	Black	1	1
19447055	Venosan Legline 20 Panythose	VL333NU	Extra Small	20mmHG	Ankle 20-21cm	Closed Toe Sheer	Nude	1	1
19447065	Venosan Legline 20 Panythose	VL335NU	Medium	20mmHG	Ankle 24-25cm	Closed Toe Sheer	Nude	1	1
19447140	Venosan Legline 30 Panythose	VL435NU	Medium	30mmHG	Ankle 24-25cm	Closed Toe Sheer	Nude	1	1
19447105	Venosan Legline 20 Panythose	VL333SA	Extra Small	20mmHG	Ankle 20-21cm	Closed Toe Sheer	Sahara	1	1
19447195	Venosan Legline 30 Panythose	VL436SA	Large	30mmHG	Ankle 26-27cm	Closed Toe Sheer	Sahara	1	1
VENOSAN 4001									
19447255	Venosan 4001 Pantyhose	SG41221	Small	18-21mmHG	Ankle 20-23cm	Closed Toe	Mexico	1	1
19447300	Venosan 4001 Pantyhose	SN41222	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19447310	Venosan 4001 Pantyhose	SN41224	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
VENOSAN 4002									
19447355	Venosan 4002 Pantyhose	SN42221	Small	23-32mmHG	Ankle 20-23cm	Closed Toe	Black	1	1
19447360	Venosan 4002 Pantyhose	SN42222	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19447365	Venosan 4002 Pantyhose	SN42223	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19447340	Venosan 4002 Pantyhose	SM42222	Medium	23-32mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19447350	Venosan 4002 Pantyhose	SM42224	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19447325	Venosan 4002 Pantyhose	SG42223	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Mexico	1	1
VENOSAN 6001									
19447460	Venosan 6002 Pantyhose	CG62212	Medium	23-32mmHG	Ankle 23-25cm	Closed Toe	Beige	1	1
19447390	Venosan 6001 Pantyhose	CG61214	Extra Large	18-21mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
VENOSAN 6002									
19447485	Venosan 6002 Pantyhose	CG62233	Large	23-32mmHG	Ankle 26-28cm	Closed Toe	Beige	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
SUPPORT SOCKS									
HEALTH SUPPORT SYSTEM									
19008120	Health Support System Sock	HSML 35-40 A R Beige	Size A	Class 3	Ankle 18-20cm	Closed Toe	Beige	1	1
19008070	Health Support System Sock	HSML 30-35 A R Black	Size A	Class 2	Ankle 18-20cm	Closed Toe	Black	1	1
19008085	Health Support System Sock	HSML 30-35 B S Black	Size B	Class 2	Ankle 20-22.5cm	Closed Toe	Black	1	1
19008095	Health Support System Sock	HSML 30-35 C S Black	Size C	Class 2	Ankle 22.5-25.5cm	Closed Toe	Black	1	1
THIGH LENGTH									
JOBST HOSIERY ULTRASHEER									
19360955	Jobst Hosiery Ultrasheer	73636-03	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360500	Jobst Hosiery Ultrasheer	75169-04	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19360505	Jobst Hosiery Ultrasheer	75170-04	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360905	Jobst Hosiery Ultrasheer	75206-04	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19361305	Jobst Hosiery Ultrasheer	75246-04	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360510	Jobst Hosiery Ultrasheer	75171-04	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360910	Jobst Hosiery Ultrasheer	75207-04	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Black	1	1
19360540	Jobst Hosiery Ultrasheer	75169-22	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19360920	Jobst Hosiery Ultrasheer	75205-23	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19360545	Jobst Hosiery Ultrasheer	75170-22	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360925	Jobst Hosiery Ultrasheer	75206-23	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19361325	Jobst Hosiery Ultrasheer	75246-20	Medium	30-40mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360550	Jobst Hosiery Ultrasheer	75171-22	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19360930	Jobst Hosiery Ultrasheer	75207-23	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19361330	Jobst Hosiery Ultrasheer	75247-20	Large	30-40mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19360555	Jobst Hosiery Ultrasheer	75172-22	Extra Large	15-20mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
19360935	Jobst Hosiery Ultrasheer	75208-23	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Natural	1	1
JOBST RELIEF HOSIERY									
19361530	Jobst Relief Hosiery	78047-09	Medium	30-40mmHG	Ankle 21-25cm	Close Toe (Pair)	Beige	1	1
19361545	Jobst Relief Hosiery	78049-08	Extra Large	30-40mmHG	Ankle 29-33cm	Close Toe (Pair)	Beige	1	1
19362200	Jobst Relief Hosiery	78042-09	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Beige	1	1
19362205	Jobst Relief Hosiery	78043-09	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Beige	1	1
19362210	Jobst Relief Hosiery	78044-09	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Beige	1	1
19362220	Jobst Relief Hosiery	78048-09	Large	30-40mmHG	Ankle 21-25cm	Closed Toe	Beige	1	1
19362215	Jobst Relief Hosiery	78045-09	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Beige	1	1
19362300	Jobst Relief Hosiery	78042-08	Small	20-30mmHG	Ankle 18-21cm	Open Toe	Beige	1	1
19362305	Jobst Relief Hosiery	78043-08	Medium	20-30mmHG	Ankle 21-25cm	Open Toe (Pair)	Beige	1	1
19362355	Jobst Relief Hosiery	78047-08	Medium	30-40mmHG	Ankle 21-25cm	Open Toe	Beige	1	1
19362310	Jobst Relief Hosiery	78044-08	Large	20-30mmHG	Ankle 25-29cm	Open Toe	Beige	1	1
19362315	Jobst Relief Hosiery	78045-08	Extra Large	20-30mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
VENOSAN 4001									
19443405	Venosan 4001 Stocking Agh	SN41603	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19443335	Venosan 4001 Stocking Agh	SM41622	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19443360	Venosan 4001 Stocking Agh	SM41602	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
19443309	Venosan 4001 Stocking Agh	SG41602	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Mexico	1	1
19443420	Venosan 4001 Stocking Agh	SN41612	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19443385	Venosan 4001 Stocking Agh	SM41613	Large	18-21mmHG	Ankle 26-29cm	Open Toe	Marokko	1	1
19443331	Venosan 4001 Stocking Agh	SG41614	Extra Large	18-21mmHG	Ankle 29-31cm	Open Toe	Mexico	1	1
VENOSAN 4002									
19443179	Venosan 4002 Stocking Agh	SN42602	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19443200	Venosan 4002 Stocking Agh	SN42622	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Black	1	1
19443180	Venosan 4002 Stocking Agh	SN42603	Large	25-32mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
19443181	Venosan 4002 Stocking Agh	SN42604	Extra Large	25-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19443170	Venosan 4002 Stocking Agh	SM42602	Medium	25-32mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
19443171	Venosan 4002 Stocking Agh	SM42603	Large	25-32mmHG	Ankle 26-29cm	Closed Toe	Marokko	1	1
19443172	Venosan 4002 Stocking Agh	SM42604	Extra Large	25-32mmHG	Ankle 29-31cm	Closed Toe	Marokko	1	1
19443308	Venosan 4002 Stocking Agh	SG41601	Small	18-21mmHG	Ankle 19-23cm	Closed Toe	Mexico	1	1
19443134	Venosan 4002 Stocking Agh	SG42601	Small Top	25-32mmHG	Ankle 19-23cm	Closed Toe	Mexico	1	1
19443183	Venosan 4002 Stocking Agh	SN42612	Medium	25-32mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19443184	Venosan 4002 Stocking Agh	SN42613	Large	25-32mmHG	Ankle 26-29cm	Open Toe	Black	1	1
19443125	Venosan 4002 Stocking Agh	SG42611	Small	25-32mmHG	Ankle 19-23cm	Open Toe	Mexico	1	1
19443305	Venosan 4002 Stocking Agh	SG41612	Medium	18-21mmHG	Ankle 23-26cm	Open Toe	Mexico	1	1
VENOSAN 5001									
19447600	Venosan 5001 Stocking Agh	SM45602	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Marokko	1	1
VENOSAN LEGLINE									
19445610	Venosan Legline 20 Agh stocking	VL325BL	Medium	20mmHG	Ankle 24-26cm	Closed Toe Sheer	Black	1	1
19445615	Venosan Legline 20 Agh stocking	VL325NU	Medium	20mmHG	Ankle 24-26cm	Closed Toe Sheer	Nude	1	1
VENOSAN 6001									
19445725	Venosan 6001 Stocking Agh	CG61621	Small	18-21mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19445710	Venosan 6001 Stocking Agh	CG61612	Medium	18-21mmHG	Ankle 23-26cm	Closed Toe	Beige	1	1
19445715	Venosan 6001 Stocking Agh	CG61613	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19445755	Venosan 6001 Stocking Agh	CG61633	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19445835	Venosan 6001 Stocking Agh	CN61633	Large	18-21mmHG	Ankle 26-29cm	Closed Toe	Black	1	1
VENOSAN 6002									
19445870	Venosan 6002 Stocking Agh	CG61613	Large	23-32mmHG	Ankle 26-29cm	Closed Toe	Beige	1	1
19445885	Venosan 6002 Stocking Agh	CG62621	Small	23-32mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19445845	Venosan 6002 Stocking Agh	CG62601	Small	23-32mmHG	Ankle 19-23cm	Open Toe	Beige	1	1
19445890	Venosan 6002 Stocking Agh	CG62622	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Beige	1	1
19445970	Venosan 6002 Stocking Agh	CN62622	Medium	23-32mmHG	Ankle 23-26cm	Open Toe	Black	1	1
19445895	Venosan 6002 Stocking Agh	CG62623	Large	23-32mmHG	Ankle 26-29cm	Open Toe	Beige	1	1
SIGVARIS STOCK									
19001090	Sigvaris Stock Cotton	50449	Large	23-32mmHG	Thigh A-G	Open Toe	Nature	1	1
THIGH LENGTH WITH BELT									
VENOSAN 4002									
19443290	Venosan 4002 Stocking Agh	SG42171	Small	23-32mmHG	Ankle 19-23cm	Open Toe	Mexico	1	1
WAIST LENGTH									
JOBST HOSIERY									
19361000	Jobst Hosiery Ultrasheer	75209-04	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Black	1	1
19361005	Jobst Hosiery Ultrasheer	75210-04	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19360605	Jobst Hosiery Ultrasheer	75174-04	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Black	1	1
19361365	Jobst Hosiery Ultrasheer	75212-04	Extra Large	20-30mmHG	Ankle 29-33cm	Closed Toe	Black	1	1
19360640	Jobst Hosiery Ultrasheer	75173-12	Small	15-20mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19361020	Jobst Hosiery Ultrasheer	75209-13	Small	20-30mmHG	Ankle 18-21cm	Closed Toe	Natural	1	1
19361025	Jobst Hosiery Ultrasheer	75210-13	Medium	20-30mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360645	Jobst Hosiery Ultrasheer	75174-12	Medium	15-20mmHG	Ankle 21-25cm	Closed Toe	Natural	1	1
19360650	Jobst Hosiery Ultrasheer	75175-12	Large	15-20mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
19361030	Jobst Hosiery Ultrasheer	75211-13	Large	20-30mmHG	Ankle 25-29cm	Closed Toe	Natural	1	1
JOBST RELIEF									
19361555	Jobst Relief Hosiery	114657A	Medium	20-30mmHG	Ankle 21-25cm	Close Toe (Pair)	Beige	1	1
19361560	Jobst Relief Hosiery	114658A	Large	20-30mmHG	Ankle 25-29cm	Close Toe (Pair)	Beige	1	1
19361575	Jobst Relief Hosiery	114662A	Large	20-30mmHG	Ankle 25-29cm	Open Toe (Pair)	Beige	1	1
19362390	Jobst Relief Hosiery	114663A	Extra Large	20-30mmHG	Ankle 29-33cm	Open Toe	Beige	1	1
VENOSAN 4002									
19443190	Venosan 4002 Stocking Ad	SG42151	Small	25-32mmHG	Ankle 20-23cm	Open Toe	Mexico	1	1
19443295	Venosan 4002 Stocking Agh	SG42172	Medium	23-32mmHG	Ankle 23-26cm	Open Toe Right	Mexico	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
COMPRESSION GARMENTS									
JOBST ELEVAREX									
19000095	Jobst Elevarex	79558-021		20mmHG				1	1
19000145	Jobst Elevarex	79558-016		20mmHG	2.5cm			1	1
19000090	Jobst Elevarex Silicone - IV	79558-015		20mmHG				1	1
JOBST ELVAREX									
19000200	Jobst Elvarex Silicone	79558-022		20mmHG	2.5cm			1	1
19000205	Jobst Elvarex Silicone - VI Plus	79558-023		20mmHG	2.5cm			1	1
19000210	Jobst Elvarex Soft	79559-023		20-30mmHG	2.5cm			1	1
19000215	Jobst Elvarex 3 Soft Ad	73579-023		30-40mmHG	2.5cm			1	1
CUSTOM MADE PRODUCT									
VENOSAN 5001									
19440015	Venosan 5001 Stocking Lymphoedema	SG45005	Custom Made	18-21mmHG			Mexico	1	1
VENOSAN 7002									
19440000	Venosan 7002 Stocking Lymphoedema	SG72505	Custom Made	23-32mmHG			Beige	1	1
19440020	Venosan 7002 Stocking Lymphoedema	SG46005	Custom Made	23-32mmHG			Mexico	1	1
VENOSAN 7003									
19440025	Venosan 7003 Stocking Lymphoedema	SG47005	Custom Made	34-46mmHG			Mexico	1	1
19440005	Venosan 7003 Stocking Lymphoedema	SG73505	Custom Made	34-46mmHG			Beige	1	1
19440010	Venosan 7003 Stocking Lymphoedema	SG73506	Custom Made	34-46mmHG			Beige	1	1
BELOW KNEE									
VENOSAN 7002									
19448055	Venosan 7002 Ad Seacell	SG72031	Small	23-32mmHG	Ankle 20-23cm	Closed Toe	Beige	1	1
19448060	Venosan 7002 Ad Seacell	SG72032	Medium	23-32mmHG	Ankle 23-25cm	Closed Toe	Beige	1	1
19448025	Venosan 7002 Ad Seacell	SG72013	Large	23-32mmHG	Ankle 25-29cm	Closed Toe	Beige	1	1
19448070	Venosan 7002 Ad Seacell	SG72034	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Beige	1	1
19448110	Venosan 7002 Ad Seacell	SN72014	Extra Large	23-32mmHG	Ankle 29-31cm	Closed Toe	Black	1	1
19448120	Venosan 7002 Ad Seacell	SN72022	Medium	23-32mmHG	Ankle 23-25cm	Open Toe	Black	1	1
19448045	Venosan 7002 Ad Seacell	SG72023	Large	23-32mmHG	Ankle 29-31cm	Open Toe	Beige	1	1
VENOSAN 7003									
19448205	Venosan 7003 Ad Seacell	SG73002	Medium	34-46mmHG	Ankle 23-25cm	Open Toe	Beige	1	1
19448320	Venosan 7003 Ad Seacell	SN73021	Small	34-46mmHG	Ankle 20-23cm	Open Toe	Black	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	MEASUREMENT	COLOUR	PACKET QTY	CARTON QTY
DONNING GLOVES							
JOBST							
19305500	Jobst Donning Gloves (Contains Latex)	79799-04	Medium			1	1
VENOSAN							
19444030	Venosan Donning Gloves	VC003	Small/Medium		White	1	10
19444035	Venosan Gloves For Doffing & Donning Nylon	VC004	Large/Extra Large		White	1	10
ROLL ON - FIXATION							
JOBST							
19364100	Jobst It Stays Roll On Body Fixation	112013	59ml			1	1
VENOSAN							
19444050	Venosan Roll On Stocking Adhesive	KC67000				1	1
STOCKING AID							
STOCKING ASSIST							
19004000	Stocking Assist Applicator	SAS		Up to 40cm Circumference	Blue	1	1
19004005	Stocking Assist	SAW		40cm		1	1
EZY AS STOCKING APPLICATOR							
19310200	Ezy As Stocking Applicator	SM259	Small	Calf 28-35cm	Orange	1	1
19310205	Ezy As Stocking Applicator	SM260	Medium	Calf 32-40cm	Yellow	1	1
19310210	Ezy As Stocking Applicator	SM261	Large	Calf 38-47cm	Blue	1	1
19310250	Ezy As Stocking Applicator Handle Attachement	SM262		Fits all Ezy As Sizes		1	1
EUREKA							
19444000	Eureka On Stocking Applicator Slide On Material	EUREKA-ON			Blue	1	1
SOCKS JET							
19444020	Socks Jet Stock Applicator W/Doffing Bars	AV002			Blue	1	1
VENOTRAIN							
19500100	Venotrain Glider Stocking Aid	26101004500000			Blue	1	1
STOCKING APPLICATOR							
EUREKA COMBI							
19447605	Eureka Combi On And Off Stocking Applicator Unisex Slide On	EUREKA-COMBI				1	1
STOCKING DONNER							
SIGVARIS DOFF N DONNER							
19000870	Sigvaris Doff N Donner Set	65246				1	1
19000871	Sigvaris Doff N Donner Cone	583D200/19000871				1	1
19000872	Sigvaris Doff N Donner Bubble	583D100				1	1
HADDENHAM GUS							
19310150	Haddenham Gus Stocking Donner	GUS	Large			1	1
JOBST							
19364000	Jobst Stocking Donner	110913				1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

STOCK CODE	PRODUCT DESCRIPTION	SUPP. PART	SIZE	COMPRESSION	MEASUREMENT	STYLE	COLOUR	PACKET QTY	CARTON QTY
ANKLE WITH ZIP									
JOBST ULCERCARE STARTER KIT									
19362700	Jobst Ulcercare Starter Kit	73631-61	Small	40mmHG	Ankle 18-21cm	Zip on Right		1	1
19362710	Jobst Ulcercare Starter Kit	73631-62	Medium	40mmHG	Ankle 21-24cm	Zip on Right		1	1
19362720	Jobst Ulcercare Starter Kit	73631-63	Large	40mmHG	Ankle 24-27cm	Zip on Right		1	1
19362730	Jobst Ulcercare Starter Kit	73631-64	Extra Large	40mmHG	Ankle 27-30cm	Zip on Right		1	1
19362740	Jobst Ulcercare Starter Kit	73631-65	XX Large	40mmHG	Ankle 30 - 33cm	Zip on Right		1	1
19362705	Jobst Ulcercare Starter Kit	73631-21	Small	40mmHG	Ankle 18-21cm	Zip on Left		1	1
19362715	Jobst Ulcercare Starter Kit	73631-22	Medium	40mmHG	Ankle 21-24cm	Zip on Left		1	1
19362725	Jobst Ulcercare Starter Kit	73631-23	Large	40mmHG	Ankle 24-27cm	Zip on Left		1	1
19362735	Jobst Ulcercare Starter Kit	73631-24	Extra Large	40mmHG	Ankle 27-30cm	Zip on Left		1	1
19362745	Jobst Ulcercare Starter Kit	73631-25	XX Large	40mmHG	Ankle 30-33cm	Zip on Left		1	1
LINERS									
JOBST ULCERCARE REPLACEMENT LINERS									
19362750	Jobst Ulcercare Replacement Liners	73632-21	Small	17mmHG	Ankle 18-21cm			3	1
19362755	Jobst Ulcercare Replacement Liners	73632-22	Medium	17mmHG	Ankle 21-24cm			3	1
19362760	Jobst Ulcercare Replacement Liners	73632-23	Large	17mmHG	Ankle 24-27cm			3	1
19362765	Jobst Ulcercare Replacement Liners	73632-24	Extra Large	17mmHG	Ankle 27-30cm			3	1
19362770	Jobst Ulcercare Replacement Liners	73632-25	XX Large	17mmHG	Ankle 30-33cm			3	1
19362775	Jobst Ulcercare Replacement Liners	73632-26	XXX Large	17mmHG	Ankle 33-36cm			3	1
19362780	Jobst Ulcercare Replacement Liners	73632-27	XXXX Large	17mmHG	Ankle 36-39cm			3	1
BELOW KNEE									
VENOSAN 8003 ULCERFIT KIT AD									
19447850	Venosan 8003 Ulcerfit Kit Ad	SG82001	Small	34-46mmHG	Ankle 20-23cm	Open Toe		1	1
19447855	Venosan 8003 Ulcerfit Kit Ad	SG82002	Medium	34-46mmHG	Ankle 23-26cm	Open Toe		1	1
19447885	Venosan 8003 Ulcerfit Kit Ad	SG82022	Medium	34-46mmHG	Ankle 23-26cm	Open Toe		1	1
19447890	Venosan 8003 Ulcerfit Kit Ad	SG82023	Large	34-46mmHG	Ankle 26-29cm	Open Toe		1	1
19447860	Venosan 8003 Ulcerfit Kit Ad	SG82003	Large	34-46mmHG	Ankle 26-29cm	Open Toe		1	1
19447930	Venosan 8003 Ulcerfit Kit Ad	SG83003	Large	34-46mmHG	Ankle 26-29cm	Open Toe		1	1
19447865	Venosan 8003 Ulcerfit Kit Ad	SG82004	Extra Large	34-46mmHG	Ankle 29-31cm	Open Toe		1	1
19447895	Venosan 8003 Ulcerfit Kit Ad	SG82024	Extra Large	34-46mmHG	Ankle 29-31cm	Open Toe		1	1
19447935	Venosan 8003 Ulcerfit Kit Ad	SG83004	Extra Large	34-46mmHG	Ankle 29-31cm	Open Toe		1	1
VENOSAN 8003 ULCERFIT AD LINER									
19447980	Venosan 8003 Ulcerfit Ad Liner	SG80003	Large	10mmHG	Ankle 26-29cm	Closed Toe	White	1	1
VENOSAN ULCERFIT AD LINER									
19447970	Venosan Ulcerfit Ad Liner	SG80001	Small	10mmHG	Ankle 20-23cm	Closed Toe	White	1	1
19447975	Venosan Ulcerfit Ad Liner	SG80002	Medium	10mmHG	Ankle 23-26cm	Closed Toe	White	1	1
19447985	Venosan Ulcerfit Ad Liner	SG80004	Extra Large	10mmHG	Ankle 29-31cm	Closed Toe	White	1	1
KNEE HIGH									
VENOTRAIN ULCERTEC STOCKING									
19500020	Venotrain Ulcertec Stocking	20980000020224	Extra Large	39mmHG	Ankle 30-33cm		White	1	1

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Product codes may vary across sizes. Products are sold individually unless otherwise specified. Prices are subject to change. See pages 7-9 for sizing charts.

Thank you

By shopping with Independence Australia, you're shopping for **a good cause.**

Independence Australia is a social enterprise that provides choices for people living with a disability or personal need, enabling them to regain and retain their independence within a supportive community.

The income we make is invested back into the vital services we provide to help people with disabilities. Your purchases help us to enhance people's lives.

We're also a registered NDIS service provider, with experts who can coordinate and develop full support packages. We provide participants with a wide range of supports and services, continence, wound care products and mobility equipment.



Shop for a
good cause

store.independenceaustralia.com
customerservice@iagroup.org.au
1300 788 855

Independence
AUSTRALIA
Your way.